

Certificate of Analysis List

For request number 677808

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2659842	4311	N/A	N/A	StablCal [®] Standard, 1.0 NTU
2660142	4288	N/A	N/A	StablCal [®] Standard, 20 NTU
2660242	4269	N/A	N/A	StablCal [®] Standard, 100 NTU
2660542	4259	N/A	N/A	StablCal [®] Standard, 800 NTU

Total Enclosures: 4

*Certificate of Analysis*

Page 1

COMMODITY: **StablCal® Standard, 1.0 NTU**COMMODITY NUMBER: **26598-42**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A4311****11/5/2014****11/10/2014**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	0.95 to 1.05 NTU	1.020 NTU

The expiration date is Nov 2016

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

*Certificate of Analysis*

Page 1

COMMODITY: **StablCal® Standard, 20 NTU**COMMODITY NUMBER: **26601-42**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A4288****10/15/2014****10/20/2014**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	19 to 21 NTU	20.0 NTU

The expiration date is Oct 2016

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist



An ISO 9001 Certified Company

P.O.Box 389
Loveland, CO 80539
(970) 669-3050*Certificate of Analysis*

Page 1

COMMODITY: **StablCal® Standard, 100 NTU**COMMODITY NUMBER: **26602-42**

MANUFACTURE DATE:

LOT NUMBER: **A4269****10/1/2014**

DATE OF ANALYSIS:

10/4/2014

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	95 to 105 NTU	100.7 NTU

The expiration date is Sep 2016

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

*Certificate of Analysis*

Page 1

COMMODITY: **StablCal® Standard, 800 NTU**COMMODITY NUMBER: **26605-42**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A4259****9/18/2014****9/22/2014**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	760 to 840 NTU	808.5 NTU

The expiration date is Sep 2016

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist