

*Certificate of Analysis*

Page 1

COMMODITY: **Formazin Turbidity Standard 4000 NTU**  
COMMODITY NUMBER: **246149**      MANUFACTURE DATE:  
LOT NUMBER: **A9022**                      **1/28/2019**

DATE OF ANALYSIS:  
**1/28/2019**

---

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity of 40 NTU formazin dilution of this lot compared to lab standard and measured on a Hach 2100 AN Turbidimeter.	39.2 to 40.8 NTU	40.70 NTU
Turbidity of 100 NTU formazin dilution of this lot compared to lab standard and measured on a Hach 2100 AN Turbidimeter.	98 to 102 NTU	102.0 NTU
Turbidity of 400 NTU formazin dilution of this lot compared to lab standard and measured on a Hach 2100 AN Turbidimeter.	392 to 408 NTU	401.0 NTU
Turbidity of 1000 NTU formazin dilution of this lot compared to lab standard and measured on a Hach 2100 AN Turbidimeter.	980 to 1020 NTU	1017.0 NTU
Turbidity of the standard undiluted.	3920 to 4080 NTU	4030.0 NTU

The expiration date is Jan 2021

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal®



*Certificate of Analysis*

Page 2

COMMODITY: **Formazin Turbidity Standard 4000 NTU**

COMMODITY NUMBER: **246149**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A9022**

**1/28/2019**

**1/28/2019**

---

***TEST***

***SPECIFICATIONS***

***RESULTS***

as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

A handwritten signature in cursive script that reads "Scott Als".

Certified by \_\_\_\_\_

Scott Als  
Analytical Services Chemist