

Certificate of Analysis

300 Technology Drive Christiansburg, VA 24073 USA inorganicventures.com

P: 800-669-6799/540-585-3030 F: 540-585-3012 info@inorganicventures.com

ACCREDITATION / REGISTRATION 1.0

INORGANIC VENTURES is accredited to ISO 17034. "General Requirements for the Competence of Reference Material Producers" and ISO/IEC 17025, "General Requirements for the Competence of Testing and Calibration Laboratories". Inorganic Ventures is also an ISO 9001 registered manufacturer (QSR Certificate Number QSR-1034).



PRODUCT DESCRIPTION 2.0

Product Code:

Single Analyte Ion Chromatography Solution

Catalog Number:

ICF1

Lot Number:

S2-F704242

Matrix:

H20

Value / Analyte(s):

1 000 µg/mL ea:

Fluoride

Starting Material:

Sodium fluoride

Starting Material Lot#:

1495, 2064, 2199

Starting Material Purity: 99.0000%

CERTIFIED VALUES AND UNCERTAINTIES 3.0

Certified Value:

 $1001 \pm 3 \, \mu g/mL$

Density:

1.001 g/mL (measured at 25 ± 4 °C)

Assay Information:

Assay Method #1

1001 ± 1 µg/mL

IC Assay NIST SRM 3183 Lot Number: 140203

Assay Method #2

1000 ± 3 µg/mL

Calculated NIST SRM Lot Number: See Sec. 4.2

- The Calculated Value is a value calculated from the weight of a starting material that has been certified directly vs. a National Institute of Standards and Technology (NIST) SRM/RM. See Sec 4.2 for balance traceability.

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

Page 1 of 5

Α	PROVADO
Respon	savel: RENATO
Padrão:	MRC Z-F1000
Data:	31/05/2022
Validad	e: 19/05/2025

Characterization of CRM/RM by Two or More Methods

Certified Value, X_{CRM/RM}, where two or more methods of characterization are used is the weighted mean of the results:

 $X_{CRM/RM} = \Sigma(w_i) (X_i)$

X_i = mean of Assay Method i with standard uncertainty uchar i

w_i = the weighting factors for each method calculated using the inverse square of

 $w_i = (1/u_{char})^2 / (\Sigma(1/(u_{char})^2))$

CRM/RM Expanded Uncertainty (±) = $U_{CRM/RM} = k (u^2_{char} + u^2_{bb} + u^2_{lts} + u^2_{ts})^{\frac{1}{2}}$

k = coverage factor = 2

 $u_{char} = [\Sigma((w_i)^2 (u_{char})^2)]^{\frac{1}{2}}$ where u_{char} i are the errors from each characterization method

ubb = bottle to bottle homogeneity standard uncertainty

uits = long term stability standard uncertainty (storage)

uts = transport stability standard uncertainty

Characterization of CRM/RM by One Method

Certified Value, X_{CRM/RM}, where one method of characterization is used is the mean of individual results:

XCRM/RM = (Xa) (uchar a)

Xa = mean of Assay Method A with

uchar a = the standard uncertainty of characterization Method A

CRM/RM Expanded Uncertainty (±) = U_{CRM/RM} = k (u²_{char a} + u²_{bb} + u²_{lts} + u²_{ts}) 1/4

k = coverage factor = 2

 u_{char} a = the errors from characterization u_{bb} = bottle to bottle homogeneity standard uncertainty

uits = long term stability standard uncertainty (storage) uts = transport stability standard uncertainty

4.0 TRACEABILITY TO NIST

- This product is traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM/RM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRM/RM are available, the term 'in-house std.' is specified.

4.1 Thermometer Calibration

- All thermometers are NIST traceable through thermometers that are calibrated by an accredited calibration laboratory.


4.2 Balance Calibration

- All analytical balances are calibrated by an accredited calibration laboratory and procedure. The weights used for testing are annually compared to master weights and are traceable to NIST.

4.3 Glassware Calibration

- An in-house procedure is used to calibrate all Class A glassware used in the manufacturing and quality control of CRM/RMs.

5.0 CHROMATOGRAM



Metrohm Compact IC Flex 930

Metrosep A Supp 5 150 x 4mm Analytical Column:

2 mM Na2CO3 Eluent:

Guard Column: Metrosep A Supp 4/5 Guard 4mm

Eluent Flow Rate: 0.7 mL/min

Anion Self Regen

Column Temp:

45°C

Suppressor/

Cell Temp:

35°C

Chemical Suppression:

Scale X-Axis:

minutes

Cation Self Regen Suppressor/

Chemical

Scale Y-Axis:

36 µS/cm

Suppression:

Concentration:

10 μg/g

Suppressor

Current/ Chemical

500 mM H2SO4

Suppressant:

6.0 INTENDED USE

- For the calibration of analytical instruments and validation of analytical methods as appropriate.

INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL 7.0

Metrohm MSM/MCS Chemical

7.1 Storage and Handling Recommendations

N/A

APROVADO

Page 3 of 5

APROVADO	
Responsável: RENATO	
Padrão: MRCZ-FICO	
Data: 31/05/2022	_
Validade: 19/05/2025	
validade: 7. 700 (2023	-

- Store between approximately 4° 30° C while in sealed TCT bag.
- While stored in the sealed TCT bag, transpiration of this CRM/RM is negligible. After opening the sealed TCT bag transpiration of the CRM/RM will occur, resulting in a gradual increase in the analyte concentration(s). It is the responsibility of the user to account for this effect. When the bottle is weighed both before and after being placed in storage, the mass difference observed will be a measure of transpiration mass loss.
- After opening the sealed TCT bag, keep cap tightly sealed when not in use and store between 4° 24° C to minimize the effects of transpiration. Use at $25^{\circ} \pm 4^{\circ}$ C to minimize volumetric dilution error when using the reported density. Do not pipette from the container. Do not return removed aliquots to container.
- For more information, visit www.inorganicventures.com/TCT

8.0 HAZARDOUS INFORMATION

- Please refer to the Safety Data Sheet for information regarding this CRM/RM.

9.0 HOMOGENEITY

- This solution was mixed according to an in-house procedure and is guaranteed to be homogeneous. Homogeneity data indicate that the end user should take a minimum sample size of 0.2 mL to assure homogeneity.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001 Quality Management System Registration

- QSR Certificate Number QSR-1034

10.2 ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories"

- Chemical Testing - Accredited / A2LA Certificate Number 883.01

10.3 ISO 17034 "General Requirements for the Competence of Reference Material Producers"

- Reference Material Producer - Accredited / A2LA Certificate Number 883.02

Inorganic Ventures, 300 Technology Drive, Christiansburg, Va. 24073, USA; Telephone: 800.669.6799; 540.585.3030, Fax: 540.585.3012; inorganic ventures.com; info@inorganic ventures.com

11.0 CERTIFICATION, LOT EXPIRATION AND PERIOD OF VALIDITY

11.1 Certification Issue Date

May 19, 2021

- The certification is valid within the measurement uncertainty specified provided the CRM/RM is stored and handled in accordance with instructions given in Sec 7.1. This certification is nullified if instructions in Sec 7.1 are not followed or if the CRM/RM is damaged, contaminated, or otherwise modified.

11.2 Lot Expiration Date

- May 19, 2025
- The date after which this CRM/RM should not be used.
- The lot expiration date reflects the period of time that the stability of a CRM/RM can be supported by long term stability studies conducted on properly stored and handled CRM/RMs. Lot expiration is limited primarily by transpiration (loss of water from the solution) and infrequently by chemical stability.

P

Responsável: RENAMO
Padrão: MRCZ-FICCO
Data: 31/05/2022
Validade: 10/05/2025

Page 4 of 5

11.3 Period of Validity

_	Sealed	TCT	Rag	Onen	Date:		
-	Sealed	101	Day	Open	Date.	Washington and the same of the	

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS Certificate Prepared By:

Uyen Truong Supervisor, Product Documentation Mya Truong

Certificate Approved By:

Michael Booth Director, Quality Control Michael 2 Booth

Certifying Officer:

Paul Gaines Chairman / Senior Technical Director

APROVADO

Responsável: RENATO

Padrão: MRC Z - F1000

Data: 31/05/2022

Validade: 19/05/2025

⁻ This CRM/RM should not be used longer than one year (or six months in the case of a 30 mL bottle) from the date of opening the aluminized bag or after the date given in Sec. 11.2, whichever comes first. This is contingent upon the CRM/RM being stored and handled in accordance with the instructions given in Sec. 7.1.