

P.O.Box 389 Loveland, CO 80539 (970) 669-3050

Certificate of Analysis This is a Component of 2971205 lot A8032

Page 1

COMMODITY: 10 NTU Verification Standard

COMMODITY NUMBER: 2961801

MANUFACTURE DATE:

DATE OF ANALYSIS:

1/30/2018

LOT NUMBER: A8012

1/18/2018

TEST

SPECIFICATIONS

RESULTS

Turbidity

9.5 to 10.5 NTU

10.13 NTU

The expiration date is Apr 2019

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

APROVADO

Responsável: RENAMO

Padrão: 29618-01

Data: 13/08/2018

Validade: 01/04/2019

Certified by

Scott Als

Analytical Services Chemist



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

Certificate of Analysis This is a Component of 2971205 lot A8032

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD

20 NTU

DATE OF ANALYSIS:

LOT NUMBER: A8019

COMMODITY NUMBER: 2684801

MANUFACTURE DATE: 1/26/2018

2/9/2018

TEST

SPECIFICATIONS

RESULTS

Turbidity

19 to 21 NTU

20.2 NTU

The expiration date is Apr 2019

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Validade: 01

Certified by .

Scott Als

Analytical Services Chemist



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

Certificate of Analysis This is a Component of 2971205 lot A8032

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD

100 NTU

DATE OF ANALYSIS:

COMMODITY NUMBER: 2684901 LOT NUMBER: A8010 MANUFACTURE DATE: 1/15/2018

1/19/2018

TEST

SPECIFICATIONS

RESULTS

Turbidity

95 to 105 NTU

102.3 NTU

The expiration date is Apr 2019

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

APROVADO

Responsável: RENAMO

Padrão: 26849-01

Data: 13/08/2018

Validade: 01/04/2019

Certified by _

Scott Als

Analytical Services Chemist



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

Certificate of Analysis This is a Component of 2971205 lot A8032

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD 800 NTU

COMMODITY NUMBER: 2660501

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: A8015 1/23/2018

1/31/2018

TEST

SPECIFICATIONS

RESULTS

Turbidity

760 to 840 NTU

812.6 NTU

The expiration date is Apr 2019

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

APROVADO

Responsável: RENAMO

Padrão: 26605-01

Data: 13/08/2018

Validade: 01/04/2019

Certified by

Scott Als Analytical Services Chemist