# Certificate of Analysis List For request number 833045

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	d Description
2659842	4358	N/A	N/A	StablCal <sup>®</sup> <b>Standard, 1.0 NTU</b>
2660101	6028	N/A	N/A	STABLCAL® FORMAZIN STANDARD 20
<b>NTU</b> 2660401 <b>NTU</b>	6040	N/A	N/A	STABLCAL® FORMAZIN STANDARD 200
2660601	6034	N/A	N/A	STABLCAL® FORMAZIN STANDARD 1000
NTU				
246102	6040	N/A	N/A	StablCal® Formazin Standard 4000
NTU				
2584201 <b>NTU</b>	6021	N/A	N/A	STABLCAL® FORMAZIN STANDARD 7500

Total Enclosures: 6



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: StablCal® Standard, 1.0 NTU

COMMODITY NUMBER: 26598-42

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: A4358

12/24/2014

12/29/2014

**TEST** 

#### **SPECIFICATIONS**

*RESULTS* 

Turbidity

0.95 to 1.05 NTU

1.020 NTU

The expiration date is Dec 2016

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by

Scott Als

Analytical Services Chemist



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD 20 NTU

COMMODITY NUMBER: 26601-01 MANUFACTURE DATE:

2/5/2016

DATE OF ANALYSIS:

LOT NUMBER: A6028 2/1/2016

> **TEST SPECIFICATIONS RESULTS**

Turbidity 19 to 21 NTU 20.5 NTU

The expiration date is May 2017

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD 200 NTU

COMMODITY NUMBER: 26604-01 MANUFACTURE DATE: DATE OF ANALYSIS:

LOT NUMBER: **A6040** 2/9/2016 2/12/2016

TEST SPECIFICATIONS RESULTS

Turbidity 190 to 210 NTU 204.3 NTU

The expiration date is May 2017

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by

Scott Als Analytical Services Chemist

LOT NUMBER: A6034



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD 1000 NTU

COMMODITY NUMBER: 26606-01 MANUFACTURE DATE:

DATE OF ANALYSIS:

2/8/2016

TEST SPECIFICATIONS RESULTS

2/4/2016

Turbidity 950 to 1050 NTU 1012.5 NTU

The expiration date is May 2017

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by

Scott Als Analytical Services Chemist



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: StablCal® Formazin Standard 4000 NTU

COMMODITY NUMBER: 2461-02 MANUFACTURE DATE:

LOT NUMBER: A6040 2/8/2016 DATE OF ANALYSIS:

2/13/2016

**TEST SPECIFICATIONS** *RESULTS* 

Turbidity 3800 to 4200 NTU 3946.2 NTU

The expiration date is May 2017

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by

Analytical Services Chemist

Scott Als



P.O.Box 389 Loveland, CO 80539 (970) 669-3050

## Certificate of Analysis

Page 1

COMMODITY: STABLCAL® FORMAZIN STANDARD 7500 NTU

COMMODITY NUMBER: 25842-01 MANUFACTURE DATE:

LOT NUMBER: A6021 1/22/2016 1/27/2016

DATE OF ANALYSIS:

7618.7 NTU

**TEST SPECIFICATIONS RESULTS** 

Formazine Concentration 7125 to 7875 NTU

The expiration date is Apr 2017